

Calprotectin Rapid Test Kit

Lateral flow rapid test cassette for detection of human calprotectin in stool samples

This kit is an immunochromatographic assay for rapid quantitative determination of human calprotectin in feces. It is for in vitro diagnostics use only.

TEST PRINCIPLE

This calprotectin Rapid Test kit contains a fecal collection tube and a test cassette. Fecal specimen is collected in the collection tube containing sample extraction buffer and then added to the test device. When the sample is added to sample pad, it moves through the conjugate pad and mobilizes colloid gold conjugated antihuman calprotectin antibody. The mixture moves along the membrane by capillary action and reacts with another anti-human calprotectin antibody that is coated on the test region. If calprotectin is present of 50ng/g stool or greater, the result is the formation of a colored band in the test region. If there is no calprotectin in the sample, the test area will remain colorless. The sample continues to move to the control area where Goat anti-mouse IgG antibody will capture colloid gold antibody conjugate to form a pink to purple color, indicating the test is working and the result is valid

REAGENTS AND MATERIALS

•	Calprotectin Cassette	REF 30823	20
•	Sampling Tube	REF 30824	20
•	Instruction for use	REE 31039	1

- Instruction for use **REF** 31039
- Rapid Cassette Test Reader (EPC 100): inquire

PRECAUTION FOR USERS

- For in-vitro diagnostic use only. 1
- Handle all specimens as if they contain infectious agents. When the assay 2. procedure is completed, dispose of specimens carefully after autoclaving for at least one hour. Alternatively, treat with a 0.5 or 1% solution of sodium hypochlorite for one hour before disposal.
- 3. Wear protective clothing (laboratory coats and disposable gloves) when assaying samples.
- Do not eat, drink or smoke in areas where specimens and kit reagents are 4. handled.
- Avoid contact between hands and eyes or nose during specimen collection 5 and testing.

SPECIMEN COLLECTION

- Collect a random sample of feces in a clean, dry container.
- Trying not to spill the buffer, unscrew and remove the white color cap with 2. attached sampling stick from the collection tube.
- Insert the sampling stick into the feces about 15 mm a few times. 3
- If necessary, remove excess of feces from the stick by gently wiping it with 4. an absorbent tissue.
- Reinsert the stick into the tube and tighten the cap thoroughly. 5
- 6. For liquid sample, add 100 µl of stool using an appropriate pipette into the vial via the green cap site.

STORAGE OF TEST KIT

The calprotectin kit can be stored at any temperature between 2-30°C. Do not freeze. The stability of the kit under these storage conditions is 18 months. Use up the reagents as soon as possible after the kit is unpacked.

TEST PROCEDURE

- Step 1. Mixing and dissolve stool sample completely by shaking or votexing.
- Step 2. Unscrew the natural colored cap on top of the green cap. Add 2 drops to the sample well of the test cassette.
- Step 3. Incubate the test at room temperature and read the test on EPC 100

Reader between 5-10 minutes.

INTERPRETATION OF RESULTS



Negative: Only one colored band appears on the control region. No apparent band on the test region.

Positive: In addition to a pink colored control band, a distinct pink colored band will also appear in the test region.

Invalid: A total absence of color in both regions is an indication of procedure error and/or that test reagent deterioration has occurred.

LIMITATIONS OF THE ASSAY

As with all diagnostic, a definitive clinical diagnosis should not be based on the result of a single test, but only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE

Sensitivity:

REF KT-112

The analytical sensitivity of the calprotectin test is 50ng/g stool. Sample containing less than 25 ng/g stool calprotectin would not show a test line.

Specificity:

The test is specific to human calprotectin. Samples containing the following substances both positive and negative controls with no effect on the test results.

	SUBSTANCES	CONCENTRATIONS
•	Bovine calprotectin	1 mg/ml
•	Goat calprotectin	1 mg/ml
•	Horse calprotectin	1 mg/ml
•	Pork calprotectin	1 mg/ml
•	Rabbit calprotectin	1 mg/ml

Interference testing:

The following substances were added to calprotectin free, 50 ng/g and 250 ng/g controls. No interference was found with any of the substances at the following concentrations.

•	Acetaminophen	20 mg/dl
•	Acetylsalicyclic Acid	20 mg/dl
•	Ampicllin	40 mg/d1
•	Ascorbic Acid	40 mg/d1
•	Atropine	40 mg/d1
•	Caffeine	40 mg/d1
•	Gentisic acid	40 mg/d1
•	Glucose	2000 mg/dl
•	Human Albumin	2000mg/d1

REFERENCES

- 1. Damms, A., Bischoff., S.C. Validation and clinical significance of a new Calprotectin rapid test for the diagnosis of gastrointestinal diseases. Int. J. colorectal Dis. 23:985-992, 2008.
- 2. Sidler, Marc, A., Leach, S. T., and Day, A. S. Fecal SI00A12 and fecal Calprotectin as noninvasive markers for inflammatory bowel disease in children. Inflamm. Bowel Dis. 14:359-366, 2008
- 3. Voganatsi A., Panutich A. Miyasaki, K.T., and Murthy R.K. Mechanism of extracellular release of human neutrophil Calprotectin complex. J. Leuk. Biol. 70:130-134.
- 4. Roseth A.G., Aadland E., Jahnsen J., et al. Assessment of disease activity in ulcerative colitis by fecal Calprotectin, a novel granulocyte marker protein. Digestion, 58:176-180, 1997.
- 5. Silberer, H., Kuppers, B., Mickisch O., et al. Fecal leukocyte protein in inflammatory bowel disease and irritable bowel syndrome. Clin. Lab 51:117-126, 2005.
- 6. Bunn S.K., Bisset W.M., Main M.J., et al. Fecal Calprotectin: validation as a noninvasive measure of bowel inflammation I childhood inflammatory bowel disease, J. Pediatr. Gastroenterol Nutr. 33:14-22, 2001.
- 7. Roseth A.G., Kristinsson J., Fagerrhol M.K. et al. Faecal Calprotectin; a novel test for the diagnosis of colorectal cancer? Scand. J. Gastroenerol. 28:1073-1076, 1993.



Epitope Diagnostics, Inc. San Diego, CA 92121, USA



MDSS GmbH Schiffgraben 41 30175 Hannover, Germany

***	Manufacturer	IVD	For <i>in vitro</i> diagnostic use only
EC REP	Authorized representative	Ĩ	Consult instructions for use
Σn	Contains sufficient for <n> tests</n>	Ť	Keep dry
REF	Catalogue Code	X	Temperature limitation
LOT	LotNumber	2	Use by
DIL	Sample diluent		6