


EDI™ Fecal C.difficile Toxin A & B ELISA Kit

Enzyme Linked Immunosorbent Assay (ELISA) for the qualitative determination of C.difficile Toxin A & B in Feces.

REF KT-836 CE IVD 

INTENDED USE

This microplate-based ELISA (enzyme linked immunosorbent assay) kit is intended for the qualitative detection of both C. difficile Toxin A and Toxin B in feces. The assay is a useful tool as an aid of detection of C. difficile infection. It is for in vitro diagnostic use only.

SUMMARY OF PHYSIOLOGY

Clostridium difficile is a common pathogen and a major cause of infectious diarrhea in hospitalized patients. While most strains produce both Toxin A and Toxin B, some only produce one or the other. These enterotoxins, which are both proinflammatory and cytotoxic, attack the mucosal lining of the intestines. If not identified and treated in a timely fashion, can result in permanent damage to the colon or colitis. This assay employs antibodies specific to both Toxin A and Toxin B.

ASSAY PRINCIPLE

This "sandwich" ELISA is designed, developed and produced for the qualitative measurement of Toxin A and Toxin B in stool specimen. The assay utilizes the microplate-based enzyme immunoassay technique by coating highly purified antibody onto the wall of microtiter wells. Controls and extracted fecal specimen are added to microtiter wells of microplate that was coated with a highly purified monoclonal anti-Toxin A and Toxin B on its wall. During the assay, the Toxin A and B will be bound to the specific antibody coated plate after an incubation period. The unbound material is washed away and another HRP-conjugated monoclonal antibodies which specifically recognizes the protein of Toxin A & B is added for further immunoreactions. After an incubation period, the immunocomplex of "Anti-Toxin A & B Capture Antibody – Toxin A & B – HRP-conjugated Anti-Toxin A & B Tracer Antibody" is formed if Toxin A or B is present in the test sample. The unbound tracer antibody and other proteins in buffer matrix are removed in the subsequent washing step. HRP conjugated tracer antibody bound to the well is then incubated with a substrate solution in a timed reaction and then measured in a spectrophotometric microplate reader. The enzymatic activity of the tracer antibody bound to Toxin A & B proteins captured on the wall of each microtiter well is directly proportional to the amount of Toxin A & B level in each test specimen.

REAGENTS: PREPARATION AND STORAGE

This test kit must be stored at 2 – 8°C upon receipt. For the expiration date of the kit refer to the label on the kit box. All components are stable until this expiration date.

1. Toxin AB Antibody Coated Microplate (30994)

Microplate coated with monoclonal anti-Toxin A and Toxin B

Qty: 1 x 96 well microplate
Storage: 2 – 8°C
Preparation: Ready to Use.

2. Anti-Toxin AB Tracer Antibody (30995)

HRP-conjugated monoclonal toxin A and B antibodies.

Qty: 1 x 12 mL
Storage: 2 – 8°C
Preparation: Ready to Use.

3. ELISA HRP Substrate (10020)

Tetramethylbenzidine (TMB) with stabilized hydrogen peroxide.

Qty: 1 x 12 mL
Storage: 2 – 8°C
Preparation: Ready to Use.

4. ELISA Stop Solution (10030)

0.5 M sulfuric acid

Qty: 1 x 12 mL
Storage: 2 – 25°C
Preparation: Ready to Use.

5. Toxin AB Negative Control (30996)

Toxin A and Toxin B negative control in a liquid bovine serum albumin based matrix with a non-azide preservative.

Qty: 1 x 1.0 mL/Vial
Storage: 2 – 8°C ; Store at <-20°C after first use.
Do not exceed 3 freeze-thaw cycles.
Preparation: Ready to Use.

6. Toxin AB Positive Control (30997)

Toxin A and Toxin B positive control in a liquid bovine serum albumin based matrix with a non-azide preservative.

Qty: 1 x 1.0 mL/Vial
Storage: 2 – 8°C ; Store at <-20°C after first use.
Do not exceed 3 freeze-thaw cycles.
Preparation: Ready to Use.

7. Concentrated Fecal Sample Extraction Buffer (30820)

Qty: 1 x 10 mL
Storage: 2 – 8°C
Preparation: 10X Concentrate. The contents must be diluted with 90 mL distilled water and mixed well before use.

8. ELISA Wash Concentrate (10010)

Surfactant in a phosphate buffered saline with non-azide preservative.

Qty: 1 x 30 mL
Storage: 2 – 25°C
Preparation: 30X Concentrate. The contents must be diluted with 870 mL distilled water and mixed well before use.

SAFETY PRECAUTIONS

The reagents must be used for in vitro diagnostic use only. Source material which contains reagents of bovine serum albumin was derived in the contiguous 48 United States. It.

was obtained only from healthy donor animals maintained under veterinary supervision and found free of contagious diseases. Wear gloves while performing this assay and handle these reagents as if they were potentially infectious. Avoid contact with reagents containing hydrogen peroxide, or sulfuric acid. Do not get in eyes, on skin, or on clothing. Do not ingest or inhale fumes. On contact, flush with copious amounts of water for at least 15 minutes. Use Good Laboratory Practices.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Precision single channel pipettes capable of delivering 10 µL, 25 µL, 50 µL, 65 µL, 100 µL, and 1000 µL.
2. Repeating dispenser suitable for delivering 100 µL.
3. Disposable pipette tips suitable for above volume dispensing.
4. Disposable 12 x 75 mm glass or plastic tubes.
5. Disposable plastic 1000 mL bottle with cap.
6. Aluminum foil.
7. Plastic microtiter well cover or polyethylene film.
8. ELISA multichannel wash bottle or automatic (semi-automatic) washing system.
9. Spectrophotometric microplate reader capable of reading absorbance at 450 nm.

SPECIMEN COLLECTION & STORAGE

Fresh fecal sample should be collected into a stool sample collection container. It is required to collect a minimum of 1-2 mL liquid stool sample or 1-2g solid sample. The collected fecal sample must be transported to the lab in a frozen condition (-20°C). If the stool sample is collected and tested the same day, it is allowed to be stored at 2-8°C. Avoid more than 3x freeze and thaw.

ASSAY PROCEDURE

1. Reagent Preparation

1. Prior to use allow all reagents to come to room temperature (20-25 °C). Reagents from different kit lot numbers should not be combined or interchanged.
2. ELISA Wash Concentrate (10010) must be diluted to working solution prior use. Please see REAGENTS section for details.
3. Concentrated Fecal Extraction Buffer must be diluted to working solution prior use. Please see REAGENTS section for details.

2. Patient Sample Preparation

Patient samples need to be diluted 1:5 with working Fecal Extraction Buffer (1x) before being measured.

1. Label a test tube (12x75 mm) or a 4 ml plastic vial.
2. With solid stool sample, take or weigh an equivalent amount (**about 250mg or 350µL for liquid feces**) with a spatula or a disposable inoculation loop. Suspend the solid/liquid stool sample with **1 mL Fecal Extraction Buffer** and mix well on a vortex mixer.
3. Centrifuge the diluted fecal sample at 3000 rpm (800- 1500 g) for 5-10 minutes. The supernatant can be directly used in the assay. As an alternative to centrifuging, let the diluted samples sit and sediment for 30 minutes and take the clear supernatant for testing.
4. Note: If the test procedure is performed on an automated ELISA system, the supernatant must be particle-free by centrifuging the sample.
5. This sample can be stored at 2-8°C up to three (3) days and below -20°C for longer storage. Avoid more than 3x freeze and thaw cycle.

3. Assay Procedure

1. Place a sufficient number of microwell strips (30994) in a holder to run Toxin AB controls and samples in duplicate
2. Test Configuration

Row	Strip 1	Strip 2	Strip 3
A	Negative Control	SAMPLE 3	SAMPLE 7
B	Negative Control	SAMPLE 3	SAMPLE 7
C	Positive Control	SAMPLE 4	SAMPLE 8
D	Positive Control	SAMPLE 4	SAMPLE 8
E	SAMPLE 1	SAMPLE 5	SAMPLE 9
F	SAMPLE 1	SAMPLE 5	SAMPLE 9
G	SAMPLE 2	SAMPLE 6	SAMPLE 10
H	SAMPLE 2	SAMPLE 6	SAMPLE 10

3. Add **100 µL** of controls (30996, 30997) and extracted patient stool samples into the designated microwell. Mix by gently tapping the plate.
4. Cover the plate with one plate sealer and aluminum foil. Incubate at **room temperature (20-25 °C) for 60 minutes**.
5. Remove the plate sealer. Aspirate the contents of each well. Wash each well **5 times** by dispensing **350 µL** of diluted wash solution (10010) into each well, and then completely aspirate the contents. Alternatively, an automated microplate washer can be used.
6. Add **100 µL** of Toxin AB Tracer Antibody (30995) to each well. Mix by gently tapping the plate.
7. Cover the plate with one plate sealer and aluminum foil. Incubate at **room temperature (20-25 °C) for 30 minutes**.
8. Remove the plate sealer. Aspirate the contents of each well. Wash each well **5 times** by dispensing **350 µL** of diluted wash solution (10010) into each well, and then completely aspirate the contents. Alternatively, an automated microplate washer can be used.
9. Add **100 µL** of ELISA HRP Substrate (10020) into each of the wells. Mix by gently tapping the plate.
10. Cover the plate with one plate sealer and aluminum foil. Incubate at **room temperature (20-25 °C) for 20 minutes**.
11. Remove the aluminum foil and plate sealer. Add **100 µL** of ELISA Stop Solution (10030) into each of the wells. Mix by gently tapping the plate.
12. Read the absorbance at **450 nm** within **10 minutes** with a microplate reader.

PROCEDURAL NOTES

1. It is recommended that all control and unknown samples be assayed in duplicate. The average absorbance reading of each duplicate should be used for data reduction and the calculation of results.
2. Keep light sensitive reagents in the original amber bottles. Store any unused antibody coated strips in the foil zip-seal bag with desiccant to protect from moisture. Exposure of the plates to humidity drastically reduces the shelf life.
3. Careful technique and use of properly calibrated pipetting devices are necessary to ensure reproducibility of the test.
4. Incubation times or temperatures other than those stated in this insert may affect the results.
5. Avoid air bubbles in the microwell as this could result in lower binding efficiency and higher CV% of duplicate readings.
6. All reagents should be mixed gently and thoroughly prior use. Avoid foaming.

INTERPRETION OF RESULTS

Visual:

1. Positive or reactive: Any sample well that is obviously more yellow than the negative control well.
2. Negative or non-reactive: Any sample well that is not obviously more yellow than the negative control well.

Note: The negative control, as well as some patient samples, may show some slight yellow color. A sample well must be obviously darker or more yellow than the negative control well, when it is interpreted as a positive result.

ELISA Reader:

1. Calculate the average absorbance for each pair of duplicate test results.
2. Calculate the cut-off
The positive cut-off and the negative cut-off are established by using following formula.

Positive Cut-Off = $1.1 \times (\text{mean extinction of negative control} + 0.10)$

Negative Cut-Off = $0.9 \times (\text{mean extinction of negative control} + 0.10)$

3. Interpret test result

- Label a test tube (12x75 mm) or a 4 ml plastic vial
- i. Positive: patient sample extinction is greater than the Positive Cut-Off
 - ii. Negative: patient sample extinction is less than the Negative Cut-Off
 - iii. Equivocal: patient sample extinction is between the Positive Cut-Off and the Negative Cut-Off.

4. Assay quality control

1. Positive control must show an average OD reading greater than 0.6.
2. Negative control should show an average OD reading less than 0.18.

Example Data and Calculation

A typical absorbance data and the resulting negative control and positive controls are represented. This absorbance must not be used in lieu of control values run with each assay.

Row	Strip 1	OD (450 nm)
A	Negative Control	0.051
B	Negative Control	0.056
C	Positive Control	0.892
D	Positive Control	0.931
E	Sample 1	0.121
F	Sample 2	0.150
G	Sample 3	0.174
H	Sample 4	0.458

The OD of negative controls and positive control meet the Internal Quality Control Standard. The Assay is valid.

Calculate the Mean OD for negative control:

$$\text{Mean}_{\text{neg.}} = (0.051 + 0.056)/2 = 0.053$$

Calculate the Positive and Negative Cut-Off Value:

$$\text{Positive Cut-Off} = 1.1 \times (0.053 + 0.10) = 0.168$$

$$\text{Negative Cut-Off} = 0.9 \times (0.053 + 0.10) = 0.138$$

$$\text{Equivocal} = 0.139 \sim 0.167$$

Interpret the Sample Result:

Sample 1 = $0.121 \leq \text{Negative COV}$ → Negative
Sample 2 = $0.150 \leq \text{Pos. COV}$; $\geq \text{Neg COV}$ → Equivocal
Sample 3 = $0.174 \geq \text{Positive COV}$ → Positive
Sample 4 = $0.458 \geq \text{Positive COV}$ → Positive

LIMITATIONS OF THE PROCEDURE

1. The results obtained with this *Toxin AB ELISA* Test Kit serve only as an aid to diagnosis and should not be interpreted as diagnostic in themselves without taking other clinical findings such as endoscopy and biopsy, etc.

QUALITY CONTROL

To assure the validity of the results each assay must include both negative and positive controls. For a valid test, the positive control must have an absorbance of at least 0.6 OD units and the negative control must be less than 0.18 OD units. We also recommend that all assays include the laboratory's own controls in addition to those provided with this kit.

CLINICAL SENSITIVITY AND SPECIFICITY

Forty stool samples were collected from normal subjects aging 23 to 64 years old. Ten fecal samples from patients with *C. difficile* infection (also positive in *C. difficile* GDH by using EpiTope Diagnostics KT-828 kit) were also collected. All fecal samples were prepared and assayed according to the procedures described in this manual. The results are summarized below.

Test Result	Normal Subjects (True Negative)	Patients (True Positive)
Negative	40	3
Positive	0	7
Total	40	10

The clinical sensitivity is **70%** ($7/10 \times 100\% = 70\%$)

The clinical diagnostic specificity is **100%** ($40/40 \times 100\% = 100\%$)

PERFORMANCE CHARACTERISTICS

Sensitivity

A positive control of *C. difficile* toxin AB at a concentration of 200 ng/ml was serially diluted with negative control down to 100 ng/ml, 50 ng/ml, 25 ng/ml, 12.5 ng/ml and 6.25 ng/ml. All these diluted samples are measured with this kit. The results showed that 12.5 ng/mL had a negative interpretation, while 25 ng/mL showed a positive test result. Therefore, the analytical sensitivity of this kit is about 25 ng/mL.

Specificity

The assay does not cross react to the following:

Helicobacter pylori, glutamate dehydrogenase (GDH) *Cryptosporidium parvum*, *Giardia lamblia*, rotavirus and adenovirus.

Reproducibility and Precision

The reproducibility of this assay is validated by measuring two samples both in a single assay of 16-replicate determinations (intra-assay) and in 10 different assays run on different dates (inter-assay).


	Intra-Assay		Inter-Assay	
Sample	1	2	1	2
Mean OD (450 nm)	0.478	0.788	0.504	0.827
Standard Deviation	0.047	0.074	0.043	0.078
CV (5%)	9.8%	9.4%	8.5%	9.5%

Interference

Two positive samples are spiked with 5% volume of interference materials to reach a final concentration shown in the table below. All samples are tested in an assay in duplicate.

	Mean OD 450 nm			
	Additive	Amt Added (mg/mL)	Sample	Sample 2
1	Test Control	-	0.352	0.235
2	Bilirubin -L	0.4	0.373	0.214
3	Bilirubin - M	2.0	0.427	0.214
4	Bilirubin - H	10.0	0.319	0.205
5	Test Control	-	0.351	0.132
6	Hb - L	0.4	0.300	0.107
7	Hb - M	2.0	0.343	0.134
8	Hb - H	10.0	0.326	0.093
9	Lipid - L	8	0.304	0.121
10	Lipid - M	40	0.331	0.105
11	Lipid -H	200	0.365	0.080

GLOSSARY OF SYMBOLS (EN 980/ISO 15223)



WARRANTY

This product is warranted to perform as described in its labeling and literature when used in accordance with all instructions. Epitope Diagnostics, Inc. DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, and in no event shall Epitope Diagnostics, Inc. be liable for consequential damages. Replacement of the product or refund of the purchase price is the exclusive remedy for the purchaser. This warranty gives you specific legal rights and you may have other rights, which vary from state to state.

REFERENCES

1. Bartlett, J.G. (2002) Clinical practice. Antibiotic-associated diarrhea. N Engl J Med. 346: 334–339
2. Kelly, C.P (1996) Immune response to *Clostridium difficile* infection. Eur J Gastroenterol Hepatol. 8: 1048–1053
3. Lyster, David M., Howard C. Krivan and Tracy D. Wilkins. *Clostridium difficile*: Its disease and toxins.
4. Borriello, S.P., FE. Barclay, P.J. Reed, A.R. Welch, J.D. Brown, and O.W. Burden. (1987) Analysis of latex agglutination test for *Clostridium difficile* toxin A(D-1) and differentiation between *Clostridium difficile* toxins A and B, and latex reactive protein. J. Clin Path. 40:573-580.

TECHNICAL ASSISTANCE AND CUSTOMER SERVICE


For technical assistance or place an order, please contact Epitope Diagnostics, Inc. at (858) 693-7877 or fax to (858) 693-7678.

This product is developed and manufactured by



Epitope Diagnostics, Inc.
7110 Carroll Road
San Diego, CA 92121, US

Please visit our website at www.epitopediagnostics.com to learn more about our products and services.



MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany

SHORT ASSAY PROTOCOL

1. Add **100 µL** of controls and extracted patient stool samples into the designated microwells.
2. Incubate at **room temperature (20-25 °C)** for **60 minutes**.
3. Wash each well five times.
4. Add **100 µL** of Toxin AB Tracer Antibody to each well.
5. Incubate at **room temperature (20-25 °C)** for **30 minutes**.
6. Wash each well five times.
7. Add **100 µL** of ELISA HRP Substrate into each of the wells.
8. Incubate at **room temperature (20-25 °C)** for **20 minutes**.
9. Add **100 µL** of ELISA Stop Solution into each of the wells.
10. Read the absorbance at **450 nm** within **10 minutes** with a microplate reader.